Hartmann Alligator Straight Ear Forceps



Like other ENT instruments these forceps are bent. There is a joint at its middle and the 2 blades open above downwards. The tip is oval pointed, serrated and slightly grooved.

Aural Hartmann forceps- One of the functions of this instrument is to clamp the arteries in the aural cavity. Once the arteries are clamped, the flow of blood is controlled. The **aural Hartmann forceps** are also useful in grasping or holding a portion in the ear. Because of the long body of the tool, the inside of the patient's aural cavity can be reached with ease. It is also used for packing the external ear, for dressing, for mopping, or for removing a foreign body or bone chips during an operation.